Monthly CeraNews

The Orthopaedic Landscape Information Newsletter, Issue January 2017

CoC can improve revision

Ceramic-on-ceramic bearings are possibly the best option for revision THA. For his study, **Wirtz** examined cumulative revision rates reported in several registers and clinical studies. In his presentation held at the 2016 AAOT Congress in Buenos Aires he pointed out that CoC bearings eliminate several revision causes, reduce the risk of re-infection and show superior results when compared to the alternatives.

READ MORE >

Noise in all bearings

In a study on noise emission from hip implants using a patient questionnaire, **Robinson et al.** found that this phenomenon is noted in hip implants with ceramic-on-ceramic as well as with metal-on-polyethylene bearings, although with different incidence. They concluded that noise apparently is an underreported phenomenon and recommend informing all patients of possible noise emission from their THA, irrespective of the bearing type.

READ MORE >

Dislocations reduced with CoC

The revision rates for late dislocation are significantly lower with ceramic-on-ceramic bearings than with other bearing couples, when the bearing diameter is larger than 28mm. **Pitto** presented his analysis of the New Zealand register data at the 2016 congress of the Japanese Hip Society. He found the best outcome with 32mm CoC bearings and hypothesised as reason the prevention of inflammatory reactions to polyethylene and metal particles leading to fluid expansion and capsule dissociation.

READ MORE >

Published by CeramTec GmbH CeramTec-Platz 1–9, 73207 Plochingen, Germany Tel.: +49 7153 611-828, Tel: +49 7153 611-950 ceranews@ceramtec.de, www.biolox.com Editorial board: • Hartmuth Kiefer • Steven Kurtz • Rocco Pitto • Robert Streicher

Is ceramic the most economical option?


When the costs of diagnosis and revision related to taper corrosion are taken into account, ceramic femoral heads can be more cost-effective than cobalt-chrome femoral heads. Regarding latest findings on these metal related issues, **Wyles et al.** calculated the overall financial burden. In their best-case scenario 0.875 % of all THA receive an adverse local tissue reaction (ALTR) work-up leading to revision in 0.11% of the patients. In the worst-case scenario, 3.5% of all THA patients receive an ALTR work-up and 1.75% are revised. They concluded that for the US health system "wholesale use (of ceramic femoral heads) in THAs may in fact provide the most economical **READ MORE** >

Medicare patients and THA bearing outcomes: Study insights

Elderly THA patients from the Medicare data base with CoP bearings show a reduced risk of dislocation, infection and mortality when compared to patients with MoP bearings. **Kurtz et al.** examined 315,784 US Medicare patients aged 65 years and older. They also found a trend towards reduced risk of revision with CoP in comparison to MoP bearings but the data did not reach statistical significance. When comparing patients with CoC and MoP bearings, there was no significant difference in risk of dislocation, revision, or mortality. However, there was a reduced risk of infection for patients with CoC bearings compared to MoP.

The Charlson comorbidity index was consistently one of the most important predictors for mortality, dislocation and revision as well as infection. Obesity was the most important risk factor for infection and the second most important factor for revision. The authors conclude that ceramic bearings are associated with lower risk of infection compared with MOP *READ MORE* >

Reasons for revision in THA



AOANJRR: Australian Orthopaedic Association National Joint Replacement Registry, Annual Report 2016 NJR: National Joint Registry for England, Wales, Northern Ireland and the Isle of Man, 13th Annual Report, 2016 Swedish Hip Arthroplasty Register, Annual Report 2016, procedures 2006–2015, p. 65 (Swedish edition) NZJR: New Zealand Joint Registry, 17 Year Report, procedures 1999–2015, p. 22

Aseptic loosening, dislocation and infection are the three most common reasons for revision in THA.

Executive Summary

Title	Ceramic Femoral Heads for All Patients? An Argument for Cost Containment in Hip Surgery			
Authors	Cody C. Wyles, Benjamin A. McArthur, Eric R. Wagner, Matthew T. Houdek, Jose H. Jimenez-Al- monte, Robert T. Trousdale			
Journal	Am J Orthopaedics 2016 September; 45(6):E362-E366			
Level of Evidence	None given.			
Summary	 Trunnionosis from modular connections of total hip arthroplasties (THA) is also an issue with metal-on-polyethylene (MoP) bearings, and can lead to increased complications such as painful THA or revisions (see also Monthly CeraNews 2_2016). The diagnosis and management of adverse local tissue reactions (ALTRs) is complex and cost-intensive. As ceramic femoral heads can mitigate this issue, a cost effectiveness model was developed by the authors. The cost estimation for an ALTR work-up was calculated following a published guideline (Kwon et al JBJS 2014). Aspects of this were imaging (MRI, ultrasonography, radiography), serum and aspiration tests, and clinical appointments and procedures including revision. The costs for the tests were taken from in-house data. The authors created two models: 1) additional cost for a ceramic femoral head and 2) cut-off value for cost effectiveness of a ceramic femoral head. Ceramic head prices were determined from 3 different practice sources for 2 different suppliers. The simulations were based on a previous finding that 7% of THA patients with MoP bearings present groin pain and that 12.5, 25 or 50% of these receive an ALTR work-up or are revised. The cost for a single ALTR work-up was \$5,085 with MRI and \$2,402. Revision with 3-day stay costs \$53,320 without perioperative medications and devices. Ceramic head extra cost was between \$500 and \$1,500. The authors concluded that their model suggests that ceramic femoral heads could be more cost-effective than cobalt-chrome alloy femoral heads. In regards to ALTRs, ceramic femoral ball heads show a superior safety profile and wholesale use in THA may in fact provide the most economical option on a societal scale. 			
Study Limitations	Use of 7% as the incidence of painful conventional THA (self-reported)			
	Based on only one clinical complication - trunnionosis			
	Cost from one finance department at one institution			
	Cost estimation only valid for the USA			
Key Messages	Trunnionosis is a serious complication, also with MoP implants.			
	Additional cost for a ceramic femoral head in the US was \$500 - \$1'500.			
	Ceramic femoral heads could be more cost-effective than CoCr femoral heads based on avoidance of the consequences of metal release.			





Executive Summary

Issue January 2017

Title	Outcomes of Ceramic Bearings After Primary Total Hip Arthroplasty in the Medicare Population
Authors	Steven M. Kurtz, Edmund Lau, Doruk Baykal, Bryan D. Springer
Journal	J Arthroplasty. Published online, 2016. DOI: 10.1016/j.arth.2016.02.054.
Level of Evidence	None given.
Summary	Kurtz et al analyze the outcome of total hip arthroplasty (THA) with ceramic-on-ceramic (CoC), ceramic-on-polyethylene (CoP) and metal-on-polyethylene (MoP) bearings for 315,784 US Medicare patients. They looked at periprosthetic joint infection (PJI), dislocation, revision, and death. Propensity scores were developed (used to treat large data sets of retrospective registry data, such as the Medicare claims administrative data) to adjust for selection bias in the choice of bearing couples. Most patients received MoP (74.7%), followed by CoP (22.3%) and CoC (3%) bearings. Patients were on average 74.3 years of age, with CoP and CoC used more often in the age group below 70 years (MoP 24.2%, CoP 40.3%, CoC 38.2% of patients). 62% were females, 94% were white, and 57% were without significant comorbidities. THA patients with CoP bearings exhibited a significantly reduced risk of dislocation (p<0.01), infection (p=0.001) and mortality (p=0.001) compared to patients with MoP bearings. Additionally a trend towards reduced revision risk with CoP in comparison to MoP bearings (p=0.095) was reported. The Charlson comorbidity index was consistently one of the most important predictors for mortality, dislocation, revision and infection. Obesity was the most important risk factor for infection and the second most important factor for revision. When comparing patients with CoC and MoP bearings, there was no significant difference in risk of dislocation, revision, or mortality. However, there was a significantly reduced risk of infection (p=0.01). The authors conclude that their study results showed no significant difference in risk of revision at 8-9 years follow up for THAs with any bearing. However, after adjusting for selection bias and various other confounding factors, ceramic bearings exhibit an association with lower risk of infection compared with MoP bearings.
Study Limitations	Analysis is limited to THAs from the Medicare records including ICD-9-CM (reporting bearing material) classification and diagnosis codes. Codes recording accuracy was not tested.
	Study with only elderly patients, > 65 years
	Several relevant factors such as differences in material (PE/XPE; type of ceramic), bearing diameter are not reported in the Medicare records.
	36 mm CoC bearings, which have shown to reduce the risk of dislocation, were only available at the end of the study period.
Key Messages	Risk of infection was lower in patients with CoP and CoC bearings compared to patients with MoP bearings.
	Risk of mortality, dislocation, and infection was lower in patients with CoP bearings compared to MoP bearings.
	No significant difference in risk of overall revision rates between different bearing surfaces.
	Charlson comorbidity index was one of the most important risk factors for mortality, dislocation, revision and infection.



Executive Summary

Title	Noise characteristics in ceramic-on-ceramic vs. metal-on-polyethylene total hip arthroplasty: a comparative study
Authors	Patrick G. Robinson, Ian Anthony, Sudeep Kumar, Bryn Jones, Andrew Stark, Roland Ingram
Journal	Hip International 2016; 26 (5):492-497. DOI: 10.5301/hipint.5000383
Level of Evidence	None given.
Summary	Robinson et al. sent a hip questionnaire (Ingram Squeaky Hip Score) and the Oxford Hip Score to 1,000 patients, of which 509 responded. Patient mean age was 63.7 years with a mean post OP follow up 33 months (6-156 months). 282 patients had ceramic-on-ceramic (CoC) and 227 metal-on-polyethylene (MoP) total hip arthroplasties (THA). In the CoC group 17% of the patients reported noise (55% clicking, 26% grinding, 19% squeaking, 17% crunching, 5% grinding, 5% popping), although the difference was not significant (p=0.054). Patients with noisy hips had an average of 5 points less in their OHS (Oxford Hip Score), however, the authors state that longer follow up is necessary to link noise to poorly functioning implants. In the CoC group 42% of noise affected patients frequently/all the time compared to 26% in the MoP group. Occasional noise was reported by 38% and 37%, respectively, and rare emission of noise in 19% and 37%, respectively. Movements causing noise were bending down and standing up, as well as taking the first few steps in both groups and squatting in the CoC group. Bending down and walking was reported to cause the loudest noise. When patients rated the effect of noise on their daily lives on a scale from 0-10 (0=no effect), the CoC group had a median score of 2 (range 0-8) and the MoP group had a median score of 1 (range 0-7). The authors found no relationship between noisy hips and BMI or femoral head size. According to the authors, noise from THA is an underreported phenomenon, which currently has been focused primarily on squeaking with CoC bearings. However, they found that it should also be considered a potential "complication" with MoP bearings. They conclude that patients should generally be forewarned of possible noise emission from their THA, irrespective of the bearing.
Study Limitations	Questionnaire based study, not validated
	Patient selection criteria not given
	61% of CoC implanted with THA components previously reported with high incidence of noise generation
	Short term study (CoC 2.5 years, MoP 3.3 years)
	Big age difference between CoC and MoP patients; age had a significant effect on noise reported
Key Messages	Noise is an underreported phenomenon of uncertain significance.
	Noise is a general issue in THA and not restricted to CoC bearings.
	Study reports squeaking even with MoP THA.
	Patients should be consented of possible noise emission irrespective of bearing surfaces.

CeramTec

ТНЕ

CERAMIC EXPERTS

Ceramic-on-Ceramic in revision hip arthroplasty

Dieter C. Wirtz Department of Orthopedics and Traumatology University Hospital Bonn Germany





BIOLOX®OPTION

AAOT 28 Nov 2016


symposium Ceramtec

Clin Orthop Relat Res (2010) 468:3070–3076 DOI 10.1007/s11999-010-1399-0

CLINICAL RESEARCH

Risk of Subsequent Revision after Primary and Revision Total Joint Arthroplasty

Kevin L. Ong PhD, Edmund Lau MS, Jeremy Suggs ScD, Steven M. Kurtz PhD, Michael T. Manley FRSA, PhD






AAOT 28 Nov 2016


revisions at 10 years postoperative

AOA NJRR 2015 – 327.151 THA 9.474 aseptic R-THA


UK-NJR 2015 – 708.311 THA 79.859 aseptic R-THA



revision reasons THA & Re-THA



Dilemma – young patients are still young at time of 1th revision






orthopädie und unfallchirurgie 🥤 bonn 5

AAOT 28 Nov 2016

symposium Ceramtec

bearings in revision arthroplasty



soft - hard

hard - hard






AAOT 28 Nov 2016



orthopädie und unfallchirurgie **bonn**


	The Journal of Arthroplasty 31 (2016) 1979–1985					
T	Contents lists available a		ý				
ELSEVIER journal	homepage: www.artl	hroplastyjourn	al.org				
Revision Arthroplasty				·			
Outcomes of Ceramic Bearing the Medicare Population	gs After Revision	n Total Hip	Arthroplasty in				
Steven M. Kurtz, PhD ^{a, *} , Edmund C. I	au, MS ^b , Doruk Bay	kal, PhD ^b , Bry	an D. Springer, MD ^c				
^a Exponent, Inc, Philadelphia, Pennsylvania ^b Exponent, Inc, Menlo Park, California ^c OrthoCarolina Hip and Knee Center, Charlotte, North Carolin		80					
	(%)	-					
	su	70			_		-
US-Medicare population	sio	co -				~	-
between 2005 and 2013. n= 31.809	evi	60				M-PE	
age > 65 yr	Ř	50					-
Note: CoC limited	bud	Ē					-
availability in US (FDA)	b m	40					_
	Usage Among Revisions	30				C-PE	
	age					C-FE	
	Us;	20					-
		10					
	Bearing			-		COC	-
	äe	o E	A	*		** - *	▲ <u>-</u>
		2004	2006	2008	2010	2012	2014
universitäts klinikum bonn				Ye	ar		orthopädie und unfallchirurgie
AAOT 28 Nov 2016		5	symposium Ceram	tec			7

Revision rate THA (only osteoarthritis)



T.R., f., 73 ys.

BPO-allergy medial + superomedial acetabular defect, stem taper not damaged







Copyright © D. Wirtz

AAOT 28 Nov 2016

symposium Ceramtec

Cementless acetabular reconstruction with C-o-C, sleeved ceramic head







orthopädie und unfallchirurgie **bonn**

AAOT 28 Nov 2016


symposium Ceramtec

10

Use of sleeved ceramic haeds

minor damage on stem taper light scratches < 0.25 mm

major damage on stem taper heavily scratched, broad truncations



BIOLOX®OPTION can be used



BIOLOX®OPTION must not be used



11

Copyright © D. Wir

AAOT 28 Nov 2016

universitäts

nikum**bonn**

symposium Ceramtec

adapter sleeves for BIOLOX[®]OPTION - literature

"large ceramic heads ... with a metal adapter sleeve have no effect on corrosion of modular taper connections"

Fretting Corrosion and Trunnion Wear — Is it also a Problem for Sleeved Ceramic Heads? Preuss R, Haeussler KL, Flohr M, Streicher, RM. Seminars in Arthroplasty 2012; 4: 251-6

"the use of the Biolox[®]OPTION system in revision hip arthroplasty demonstrates little damage to either the titanium adapter sleeve or the ceramic head"

Performance of Biolox Delta Ceramic Bearings with Titanium Adapter Sleeves in Revision Hip Arthroplasty: A Retrieval Analysis. Figgie M. Jr., Elpers, Padgett D. Abstract ORS 2015

"fretting corrosion in sleeved ceramic heads showed lower levels than observed in prior studies of tapers in CoCr femoral heads. None of the sleeves in this study had severe corrosion of the internal sleeve surface"

Fretting and Corrosion Damage in Taper Adapter Sleeves for Ceramic Heads: A Retrieval Study. MacDonald DW, Chen A, Lee GC, Klein GR, Mont MA, Kurtz SM. Sumitted to JoA August 2016





AAOT 28 Nov 2016

symposium Ceramtec

ceramic-inlays in cup revision



CombiCup R- Link



MRSC - Brehm



Copyright © D. Wirt:

AAOT 28 Nov 2016


universitäts klinikum**bonn**

symposium Ceramtec

13

K.A., f., 64ys.

2.- Re-THA, BPO+Nickel allergy





orthopädie und unfallchirurgie **bonn**




AAOT 28 Nov 2016

symposium Ceramtec

14

B.H., m., 62 ys.

extended PE + ME-granuloma superomedial + craniolateral defect




AAOT 28 Nov 2016

symposium Ceramtec



"augment-and-modular cage" MRSC [Brehm]






Copyright © D. Wirtz

AAOT 28 Nov 2016

symposium Ceramtec

16

Why ceramic in hip revision? aseptic loosening, young patient



G.R., f., 64ys.

2. Re-THA, aseptic cup loosening, multiple dislocations (head 28mm)










AAOT 28 Nov 2016

symposium Ceramtec

Why ceramic in hip revision ? dislocation






"...he reasons may be related to observed differences in the periarticular muscles (fat atrophy or not) ..."

orthopädie und unfallchirurgie



US-Medicare Register 2005 - 2013

- 31'809 Revisions-THA
- adjusted for patient-, hospital- und surgeon riskfactors; Cox-Regression
- Ce/Ce significant better
 - HR = 0,76; *p*=0,04 im Vergleich zu Me/PE





Wirtz

Copyright © D.



A.M., m., 59ys

Girdlesone-situation 6 mo after explantation MRSA, Proprioni acnes





orthopädie und unfallchirurgie **(bonn** 21

AAOT 28 Nov 2016

symposium Ceramtec

Why ceramic in hip revision? Periprosthetic joint infections

Clin	Orthop Relat Res
DOI	10.1007/s11999-016-4916-y

Clinical Orthopaedics and Related Research®

SYMPOSIUM: 2015 INTERNATIONAL HIP SOCIETY PROCEEDINGS

Periprosthetic Joint Infection in Hip Arthroplasty: Is There an Association Between Infection and Bearing Surface Type?

0.35%

Rocco P. Pitto MD, PhD, Laurent Sedel MD



universitäts

inikum**bonn**




Fig. 1 This figure shows the percent of revision with CIs for PJI within 6 months after the index procedure by bearing surface.






Fig. 2 The Kaplan-Meier survival analysis shows the proportion of revision-free THAs for PJI by bearing surface. The median observation period in this patient population (84,894 THAs) was 9 years (range, 1-15 years).

Why ceramic in hip revision? Periprosthetic joint infection



- Bozic KJ, Ong K, Lau E, Kurtz SM, Vail P, Rubash H. Risk of Complication and Revision Total Hip Arthroplasty Among Medicare Patients with Different Bearing Surfaces. CORR 2010;468;2357-2362
- Trebse R, Levasic V, Kovac S. Prostethic Joint Infections and bearings. Hip International 2014;24(5), 533
- Alijanipour P, Restrepo C, Smith L, Parvizi J, Malkani A. Periprosthetic joint infection: Could Bearing Surface Play a Role? Presentation 45th Annual Meeing Eastern Orthopaedic Association 2014
- Smith L, Alijanipour P, Restrepo C, Maltenfort M, Parvizi J, Malkani A. Periprosthetic joint infection: Could Bearing Surface play a Role? Abstract, 45th Meeting of the Eastern Orthopaedic Association 2014, 197
- 5. $12^{\text{th}}\text{Annual Report: National Joint Registry for England, Wales and Northern Ireland, 2015$

- Falcioni S, Ancarani C, Bordini B, Pichierri M, Stea S. Influence of articular coupling on septic loosening of total hip arthroplasty. Abstract EHS 2014
- Varnum C, Pedersen AB, Kjaesgaard-Andersen P, Overgaard S. Comparison of the risk of revision in cementless total hip arthroplasty with ceramic-on-ceramic and metal-on-polyethylene bearings. Acta Orthopaedica 2015;86(3)
- Graves SE, Lorimer M, Bragdon C, Muratoglu O, Malchau H. Reduced risk of revision for infection when a ceramic bearing surface is used. Abstract ISTA 2015
- 9. Pitto RP, Sedel L. Periprosthetic Joint Infection in Hip Arthroplasty: Is There an Association Between Infection and Bearing Surface Type? Clin Orthop Relat Res 2016;DOI 10.1007/s1999-016-4916-y





Streicher R, unpublished data Ceramtec

Key Engineering Materials Vol. 631 (2015) pp 448-453 © (2015) Trans Tech Publications, Switzerland doi: 10.4028/www.scientific.net/KEM.631.448 Submitted: 28.06.2014 Accepted: 22.07.2014


Bioceramic materials show reduced pathological biofilm formation

Corrado PICONI^{1,a}, Andrei C. IONESCU^{2,b}, Andrea COCHIS^{3,4,c}, Erica IASI^{4,d}, Eugenio BRAMBILLA^{2,e} and Lia RIMONDINI^{3,4,f*},

¹ Centre for New Materials and Prosthetic Technologies, Università di Tor Vergata, via Montpellier 1, 00133 Rome (RM), Italy

²Department of Medical, Surgical and Dental Sciences, Università di Milano, 20121 Milan (MI), Italy

³Consorzio Interuniversitario per la Scienza e Tecnologia dei Materiali (INSTM), Via G. Giusti 2, 50121 Firenze (FI), Italy



(B)

S. epidermidis		
Specimen	% Inhibition (vs control)	
BIOLOX®forte	38.6 ± 3.3	
Si3N4	25.9 ± 4.5	
BIOLOX [®] delta	26.8 ± 3.9	


(D)

E. coli		
Specimen	% Inhibition (vs control)	
BIOLOX®forte	62.3 ± 3.5	
Si3N4	61.5 ± 3.8	
BIOLOX [®] delta	57.5 ± 3.2	



K.M., f., 63 ys.,

ME-ME resurfacing with aseptic cup loosening 4 ys. postop revision with CE-CE + cementless stem + 36 head






orthopädie und unfallchirurgie **bonn** 25

AAOT 28 Nov 2016

symposium Ceramtec

Why ceramic in hip revision ? ALRT - "wear disease"



3-body wear-study



Porporati, Laborstudie Endolab, ISTA 2014

Copyright © D. Wirtz

unfallchirurgie 🧧 bonn

B.M., f., 55 ys.

rim fracture of the ceramic inlay with multiple ceramic particles within the soft-tissues



preop





postop





AAOT 28 Nov 2016

symposium Ceramtec

Why ceramic in hip revision? Ceramic fracture





Hip International / Vol. 12 no. 1, 2002 / pp. 37-42

Case report

Massive wear in a CoCrMo head following the fracture of an alumina head

P.A. GOZZINI', C. SCHMID², P. DALLA PRIA³

¹ Orthopaedic Department, A.O.M. Mellini, Chiari (BS) - Italy ² Dept. Eng. Materials and App. Chem. (DIMCA), University of Trieste, Trieste - Italy ³ Lima-Lto SpA, Villanova di San Daniele del Friuli, Udine - Italy



The Journal of Arthroplasty Vol. 24 No. 8 2009

Third-Generation Ceramic-on-Ceramic Bearing Surfaces in Revision Total Hip Arthroplasty

Jun-Dong Chang, MD, PhD,* Rutuj Kamdar, MS,* Je-Hyun Yoo, MD, PhD,* Mina Hur, MD, PhD,† and Sang-Soo Lee, MD, PhD*



© Wichtig Editore, 2002





indications of C-o-C in Re-THA

- revision because of aseptic loosening in young patients (<70ys)
 - => avoiding PE-wear / osteolysis
- revision because of dislocation with small heads
 "upgrade" to larger head diameter (36 / 40)
- revision because of periprosthetic joint infection
 => reducing risk of re-infection
- revision because of ALTR in case of MoM
 => stops effect and eliminates the risk of re-occurrence
- revision because of ceramic fracture
 - => best and safest option







ceramic on ceramic is an (the) option in revision hip arthroplasty





Copyright © D. Wirtz

AAOT 28 Nov 2016

symposium Ceramtec

Ceramic Bearings and Revision Risk for Late Dislocation

R.P. Pitto, MD, PhD

Department of Orthopaedic Surgery Middlemore Hospital and Manukau Surgery Centre University of Auckland, New Zealand





The University of Auckland

Acknowledgements



Prof. Laurent Sedel Dr. Chris Frampton The NZ Arthroplasty Registry

Disclosure

CORR CeramTec





Copyright © R. Pitto

Introduction

- Dislocation is a major complication after THA and constitutes a prominent reason for revision surgery
- Multiple studies have focused on risk factors for late dislocation, but bearing surface has not been extensively studied
- Currently there is conflicting evidence about bearing surfaces and dislocation rates





Introduction

The incidence of dislocation is highest in the first year after arthroplasty, and then continues at a constant rate over the life of the implant

- Early (<1yr) versus late dislocation (>1yr):
 - Early: patient and surgical factors
 - Late: biological factors







 To determine whether the bearing surface is a risk factor for late revision due to dislocation in primary THA





Methods

NZ Arthroplasty Registry:

- Since 1999
- 100,315 primary THJRs (16-year Report)
- Dislocation most common cause for revision (26.5%)





Methods

- **Exclusion criteria:**
 - Resurfacing arthroplasty
 - Revision arthroplasty
- Primary endpoint was revision for late dislocation (late defined as >1 year postOP)

Material and Methods

- 73,386 THRs fulfilling inclusion:
 - − 73,386 hips ≥1 year postOP
 - 65,387 hips <u>></u>2 years postOP
 - 42,086 hips <u>></u>5 years postOP
 - 12,967 hips <u>>10 years postOP</u>
- Mean age 68.9 years
- Mean 10-year Follow-up
- 53.2% female
- 88% OA

Material and Methods

- Surgical Approach:
 - Posterior 65.3%
 - Lateral 28.3%
 - Anterior 4.2% (other 2.2%)
- Bearing surfaces:
 - MoP 53,331
 - CoP 14,093
 - CoC 8,177
 - MoM 5,910
 - CoM 461

Results

- 3130 (4.3%) hips revised for any cause
 Rate of 0.7/100 component / years
- 836 (1.1%) revised for dislocation
 Rate of 0.19/100 component / years
- 470 (0.65%) revised for dislocation >1 year
 Rate of 0.11/100 component / years

Multivariate Analysis

- **CoC HR versus**
- CoP

(HR 2.10; 95% CI 1.12 – 3.94, p=0.021)


-MoP

(HR 1.76; 95% CI 0.94 - 3.28, p=0.075)

Adjusted for age, gender, head size, surgical approach

Results

- There were statistically significant lower rates of revision for dislocation in all age groups with >28mm CoC bearings than:
 - MoM (HR = 0.36; 95% CI 0.20 0.67, p= 0.004)
 - CoP (HR = 0.51; 95% CI 0.30 0.89, p= 0.018)
 - MoP (HR = 0.55; 95% CI 0.33 0.93, p= 0.027)




Revision for dislocation after 1 year

Revision for dislocation after 1 year

Head size >28mm and age <65 years

Head size >28mm and age ≥ 65 years



There were more revisions in CoC than MoM THAs in younger patients and smaller head size (< 65 years, 28 mm) (HR 0.29; 95% CI 0.12–0.71; p = 0.014)

Discussion

- This 10-year Registry analysis shows low rates of revision for late dislocation with CoC THRs
- Confirms Australian Registry finding regarding the increased risk of revision for late dislocation in patients younger than 65 yrs with 28mm CoC
- Confirms findings of previous paper showing low rates of late dislocation with 32mm CoC

Sexton SA et al.: CoC and risk or revision due to dislocation after THA. JBJS 91B: 1448-53, 2009

Hernigou P et al.: CoC Decreases the Long-term Risk of Dislocation. CORR 471:3875-82, 2013

Discussion

- Late dislocations may be influenced by biological factors:
 - analysis of tissue reaction to ceramics has shown small numbers of macrophages, few foreign body type giant cells and occasional lymphocytes
 - polyethylene implants promote extensive foreign body type inflammatory changes

Esposito C. et al.: Periprosthetic Tissues from CoC THAs. J Arthroplasty 2013;28:860-6

Discussion

Late dislocations may be influenced by biological factors like *Pathology of the Pseudo-Capsule:*

- MoP pseudocapsules exhibit significantly higher levels of inflammatory markers than CoC
- inflammatory reaction to polyethylene and metal wear particles results in fluid expansion and capsule dissociation

Sedel L. et al.: Prostaglandin levels in peri-THA tissues. Arch Orthop Trauma Surg 1992;111:255-8

Message to Take Home

Ceramic Bearings:

Low rates of revision for late dislocation
Best outcome with 32mm bearing surfaces





